

ELISA kits available from ADI (see details at the web site)

#0010	Human Leptin		
#200-120-AGH	Human globular Adiponectin (gAcrp30)		
#0700	Human Sex Hormone Binding Glob (SHBG)		
#0900	Human IGF-Binding Protein 1 (IGFBP1)		
#1000	Human C-Reactive Protein (CRP)		
#100-110-RSH	Human Resistin /FIZZ3		
#100-140-ADH	Human Adiponectin (Acrp30)		
#100-160-ANH	Human Angiogenin		
#100-180-APH	Human Angiopoietin-2 (Ang-2)		
#100-190-B7H	Human Bone Morphogenic Protein 7 (BMP-7)		
#1190	Human Serum Albumin	#1200	Human Albumin (Urinary)
#1750	Human IgG (total)	#1760	Human IgM
#1800	Human IgE	#1810	Human Ferritin
#1210	Human Transferrin (Tf)	#0020	Beta-2 microglobulin
#1600	Human Growth Hormone (GH)		
#0060	Human Pancreatic Colorectal cancer (CA-242)		
#1820	Human Ovarian Cancer (CA125)	#1830	Human CA153
#1840	Human Pancreatic & GI Cancer (CA199)		
#1310	Human Pancreatic Lipase		
#1400	Human Prostatic Acid Phosphatase (PAP)		
#1500	Human Prostate Specific Antigen (PSA)	#1510	free PSA (fPSA)
#0500	Human Alpha Fetoprotein (AFP)		
#0050	Human Neuron Specific Enolase (NSE)		
#0030	Human Insulin	#0040	Human C-peptide
#0100	Human Luteinizing Hormone (LH)		
#0200	Human Folicle Stimulating Hormone (FSH)		
#0300	Human Prolactin (PRL)		
#0400	Human Chorionic Gonadotropin (HCG)	#0410	HCG-free beta
#0600	Human Thyroid Stimulating Hormone (TSH)		
#1100	Human Total Thyroxine (T4)	#1110	Human Free T4 (fT4)
#1650	Human free triiodothyronine (fT3)	#1700	Human T3 (toal)
#1850	Human Cortisol	#1860	Human Progesterone
#1865	Human Pregnonolone	#1875	Human Aldosterone
#1880	Human Testosterone	#1885	Human free Testosterone
#1910	Human Androstenedione	#1920	Human Estradiol
#1925	Human Estrone	#1940	Dihydrotestosterone (DHT)
#1950	Human DHEA-sulphate (DHEA-S)		
#3400	Human serum Neopterin		
#3000	Human Rheumatoid Factors IgM (RF)		
#3100	Human anti-dsDNA		
#3200	Anti-Nuclear Antibodies (ANA)		

Alpha Diagnostic Intl. (www.4adi.com) 3000-110-RFG/180627AA **Page 7**

Instruction Manual No. M-3000-110-RFG

Rheumatoid Factor IgG (RF IgG) ELISA KIT

Cat. No. 3000-110-RFG, 96 Tests

**For Quantitative Determination of RF IgG
In Human Serum or plasma**

For In Vitro Research Use Only



6203 Woodlake Center Drive • San Antonio • Texas 78244 • USA.

Phone (210) 561-9515 • Fax (210) 561-9544

Toll Free (800) 786-5777

Email: service@4adi.com

Web Site: www.4adi.com

Draft version-Please consult the manual supplied with the kit for any lot specific change.

Rheumatoid Factor IgG (RF IgG) ELISA KIT # 3000-110-RFG

For Quantitative Determination of RF IgG In Human Serum


Kit Contents: (reagents for 96 tests)

Components	
Purified gamma globulin coated microwell strips (96 wells), Ready-to-use, P-3000111	1 Plate
RF IgG Standard A (0 U/ml) 1.5 ml, #3000112A	1 vial
RF IgG Standard B (15 U/ml) 1.5 ml, #3000112B	1 vial
RF IgG Standard C (50 U/ml) 1.5 ml, #3000112C	1 vial
RF IgG Standard D (150 U/ml) 1.5 ml, #3000112D	1 vial
RF IgG Standard E (500 U/ml) 1.5 ml, #3000112E	1 vial
RF IgG Positive control , 1.5 ml, #P3000113	1 vial
RF IgG Negative control , 1.5 ml, #N3000114	1 vial
RF IgG Sample Diluent (5X), 20 ml, #3000115	1 bottle
Anti-hIgG HRP Conjugate, 15 ml, #3000116	1 bottle
Wash buffer (50X), 20 ml, #3000110-WB dilute 1:50 with distilled water,	1 bottle
HRP Substrate Solution, 15 ml, #3000110-TMB	1 bottle
Stop solution, 15 ml, #3000110-ST	1 bottle
Complete Instruction Manual# M-3000-110-RFG	1

Intended use:

ADI's Rheumatoid Factor IgG is for Quantitative Determination of RF IgG In Human Serum or plasma. For in-vitro research use only (RUO).

Introduction



Rheumatoid factor (RF) is the autoantibody (antibody directed against an organism's own tissues) that was first found in rheumatoid arthritis. It is defined as an antibody against the Fc portion of IgG (an antibody against an antibody). RF and IgG join to form immune complexes that contribute to the disease process. Although predominantly encountered as IgM, rheumatoid factor can be of any isotype of immunoglobulins, i.e. IgA, IgG, IgM, [2] IgE, [3] IgD. The presence of IgM rheumatoid factor (RF) in the serum is the sole serological indicator included in the ACR list of criteria for the diagnosis of

RA. RFs are a subset of antiglobulins directed against the Fc region of IgG. In this definition we do not include antibodies to the IgG Fab region and pepsin agglutinators, directed against neoantigens on IgG exposed by pepsin cleavage. It is claimed that the majority of antiglobulin activity in normal serum is Fab-specific, whereas an- t globulin from RA patients is mostly Fc-specific. RFs are present in the serum of 75-80% of patients with RA at some time during the disease course. However, RFs are also found in the serum of patients with infectious and autoimmune diseases, hyperglobulinemia, B-cell lymphoproliferative disorders and in the aged population. This suggests that RF may be a finding associated with B-cell hyperactivity.

Rheumatoid factors which have been found among the IgM, IgG and IgA classes of immunoglobulins, reacting only with xenogeneic Fc are not autoantibodies and are unlikely to be of pathological significance. However, RFs can bind IgG from many species, including autologous IgG, when immobilised on surfaces. Autologous binding is of higher affinity than xenogeneic binding. The here presented test systems for the determination of rheumatoid factors uses only human Fc fragments as coated antigen. It is generally considered that high RF titers are associated with more severe disease and the presence of extra-articular features and rheumatoid nodules. This conclusion may depend on the disease duration. Serum IgM RF may precede the onset of RA by several years. A high titer of RF in non-RA individuals is associated with increased risk of developing RA. In the first 2 years of RA (early RA), serum levels of IgM, IgG and IgA RF do not correlate with disease activity. Serum IgG and IgA RF in these years are prognostic of erosive

PERFORMANCE CHARACTERISTICS

PRECISION

Intra-assay precision:

Sample	Mean (U/ml)	CV%
1	14.7	5.1
2	94.5	3.6
3	308.2	4.6

Inter-assay precision:

Sample	Mean (U/ml)	CV%
1	15.9	7.1
2	82.5	3.7
3	305	2.4

EXPECTED VALUES

It is recommended that each laboratory must determine its own negative and positive values. Samples containing less than 25 IU/ml RF IgG can be considered as RF IgG-negative; samples showing greater than 25 IU/ml concentrations can be considered as RF IgG-positive. Concentrations of higher than 75 IU/ml RF IgG usually indicate rheumatoid arthritis. RF IgG are found in about 2-10% of apparently healthy Caucasian adults and in about 50-70% of adults with classical rheumatoid arthritis.

Sensitivity

The lower detection limit for Rheumatoid Factor IgG was determined at 1.0 U/ml.

Specificity

The microplate is coated with the Fc fragment of highly purified human Immunoglobulin G. The test kit is specific for all classes of rheumatoid factors.

Calibration

The quantitative test system for Rheumatoid Factor IgG is calibrated in relative arbitrary units. The calibration is related to the 1st British Standard Preparation 64/2. This material tests positive for IgG rheumatoid factors.

Species Crossreactivity


This kit is recommended for human samples only. Its utility in other species such as mouse, rat, or monkey etc has not been tested. ADI has a separate RF ELISA kit for mouse samples (#6200).

References: Johnson PM (1976) Clin. Immunol. 6, 414; Bartfield H (1969) Ann. NY Acad. Sci. 168, 30; Karah J (1980) J Immunol. Meth. 32, 115; Moors TL (1978) Arth Rheumat. 21, 935; Seymour P (1980) Am. Soc. Clin. Pathol. 74, 776; Puente A (1988) Arth. Rheumat. 31, 1230.

WORKSHEET OF TYPICAL ASSAY

Wells	Stds (U/ml)	Mean A _{450 nm}	Calcul. Conc. (U/ml)*
A1, A2	0	0.048	
B1, B2	15	0.196	
C1, C2	50	0.464	
D1, D2	150	1.046	
E1, E2	500	1.885	

NOTE: These data are for demonstration purpose only. A complete set of negative, positive, and calibrator standards set must be run in every assay to determine sample values. Each laboratory should determine their own normal reference values.



*5_ADI_ELISA-graph


CALCULATION OF RESULTS

Calculate the mean absorbance for each duplicate. Subtract the absorbance of the sample diluents from the mean absorbance values of negative & positive controls, calibrator, and samples. The RF IgG values for the patient samples can be calculated as follows:

$$\text{IU/ml of sample} = \frac{\text{Net absorbance of test sample}}{\text{Net absorbance of the calibrator}} \times \text{IU/ml of calibrator}$$

joint disease. In established RA, high titer serum IgM RF correlates with the presence of articular disease and nodules but not with systemic disease activity. The presence of either IgG or IgA RF in patients with long-standing RA may be a good prognostic indicator of systemic manifestations. IgG and IgM RF are associated with extra-articular RA including rheumatoid vasculitis and nodules. The presence of IgM RF containing immune complexes with bound complement (C1q) is also associated with extra-articular RA.

PRINCIPLE OF THE TEST



Rheumatoid Factor IgG (RF IgG) ELISA kit is based on binding of RF IgG from serum samples to human gamma globulin immobilized on microtiter wells. After a washing step, anti-human IgG-HRP conjugate is added. After another washing step, to remove all the unbound enzyme conjugate, chromogenic substrate is added and color developed. The enzymatic reaction (color) is directly proportional to the amount of RF IgG present in the sample. Adding stopping solution terminates the reaction. Absorbance is then measured on a microtiter well ELISA reader at 450 nm. and the concentration of RF IgG in samples is calculated on the basis of the absorbance of the negative, positive, and, calibrator controls.

MATERIALS AND EQUIPMENT REQUIRED

Adjustable micropipet (20-100 µl) and multichannel pipet with disposable plastic tips. Reagent troughs, plate shaker (orbital shaker), plate washer (recommended) and ELISA plate Reader.

PRECAUTIONS

The Alpha Diagnostic International Rheumatoid Factor IgG ELISA Kit is intended for *in vitro research* use only. The reagents contain thimerosal as preservative; necessary care should be taken when disposing solutions. The Negative, Positive, and Calibrator controls have been prepared from human sera shown to be negative for HBsAg and HIV antibodies. Nevertheless, such tests are unable to prove the complete absence of viruses, therefore, sera should be handled with appropriate precautions.

Calibrators, Sample Diluent, and Antibody HRP contain bromonitrodioxane (BND: 0.05%, w/v). Stop Solution contains dilute sulfuric acid. Follow good laboratory practices, and avoid ingestion or contact of any reagent with skin, eyes or mucous membranes. All reagents may be disposed of down a drain with copious amounts of water. MSDS for TMB, sulfuric acid and BND can be requested or obtained from the ADI website: Sample Diluent and anti-Protein G-HRP contain Proclin 300 (0.05%, v/v). <http://4adi.com/objects/catalog/product/extras/ELISA-Kit-SDS-MSDS-Set-1.pdf>

SPECIMEN COLLECTION AND HANDLING

Collect blood by venipuncture, allow to clot, and separate the serum by centrifugation at room temperature. Do not heat inactivate the serum. If sera cannot be immediately assayed, these could be stored at -20°C for up to six months. Avoid repeated freezing and thawing of samples. No preservatives should be added to the serum.

REAGENTS PREPARATIONS:

Wash buffer is supplied as 50x stock. Dilute 20 ml into 980 ml de-ionized or distilled water, mix, and store at room temp for 1-2 weeks. It can be stored at 4°C for long term storage.

Sample Diluent (5X): Dilute 20 ml into 80 ml de-ionized or distilled water. Dilute serum sample 1:100 in 1x sample diluent (5 µl sample in 495 µl buffer).

STORAGE AND STABILITY

The microtiter well plate and all other reagents are stable at 2-8°C until the expiration date printed on the label. The whole kit stability is usually 6 months from the date of shipping under appropriate storage conditions.

TEST PROCEDURE (ALLOW ALL REAGENTS TO REACH ROOM TEMPERATURE BEFORE USE).

Label or mark the microtiter well strips to be used on the plate. Dilute serum samples 1:100 (5 µl of sample in a total volume of 500 µl of sample diluents). **Dilute wash buffer (1:50) with distilled water (20 ml stock in total of 1-liter). Dilute Sample Diluent (5X): Dilute 20 ml into 80 ml de-ionized or distilled water. Standards and controls are ready-to-use.**

1. Pipet **100 µl** of **ready-to-use standards**, negative & positive controls, and diluted serum samples into appropriate wells in *duplicate*. Cover the plate and incubate for **30 minutes** at **room temperature (20-28°C)**.
2. Aspirate and wash the wells **3 times** with 300 µl of diluted wash buffer. We recommend using an automated ELISA plate washer for better consistency. Failure to wash the wells properly will lead to high blank or zero values. If washing manually, plate must be tapped over paper towel between washings to ensure proper washing.
3. Add **100 µl** of antibody-enzyme conjugate into **each well**. Mix gently. Cover the plate and incubate for **15 minutes** at room temperature(**20-28°C**).
4. Aspirate and wash the wells **3 times** with 300 µl of diluted wash buffer, as above.
5. Dispense **100 ul TMB substrate per well**. Mix the plate gently for 5-10 seconds. Cover the plate and incubate for **15 minutes** at room temperature. Blue color develops into standards and positive samples.
6. Stop the reaction by adding **100 µl of stop solution to all wells** and incubate for **5 minutes** at room temperature.Mix gently. Blue color turns yellow.
7. Measure the absorbance at 450 nm using an ELISA reader.

NOTES

Read instructions carefully before the assay. Do not allow reagents to dry on the wells. Careful aspiration of the washing solution is essential for good assay precision. Since timing of the incubation steps is important to the performance of the assay, pipet the samples without interruption and it should not exceed five minutes to avoid assay drift. If more than one plate is being used in one run, it is recommended to include a set of negative & positive standards and calibrator on each plate. The unused strips should be stored in a sealed bag at 4°C. Addition of the HRP substrate solution starts a kinetic reaction, which is terminated by dispensing the stopping solution. Therefore, keep the incubation time for each well the same by adding the reagents in identical sequence. Plate readers measure absorbance vertically. Do not touch the bottom of the wells.

QUALITY CONTROL OF THE TEST

1. OD values will vary with the temperature and length of incubation.
2. The O.D. of the Standards A (reagents blank) should be <0.250 and Standard E >1.100.
3. The value of positive and negative controls should be within the range of indicated value.

Calculation of results

For Rheumatoid Factor IgG a 4-Parameter-Fit with lin-log or lin-lin coordinates for optical density and concentration is the data reduction method of choice.

Recommended Lin-Log Plot

First calculate the averaged optical densities for each calibrator well. Use lin-log graph paper and plot the averaged optical density of each calibrator versus the concentration. Draw the best fitting curve approximating the path of all calibrator points. The calibrator points may also be connected with straight line segments. The concentration of unknowns may then be estimated from the calibration curve by interpolation.

Interpretation of results

In a normal range study with serum samples from healthy blood donors the following ranges have been established with the Rheumatoid Factor IgG test:

Rheumatoid Factor IgG [U/ml]	
Normal:	< 20
Elevated:	≥ 20

Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually. It is recommended that each laboratory establishes its own normal and pathological ranges of serum Rheumatoid Factor.

LIMITATIONS OF PROCEDURE

The absence of rheumatoid factor does not rule out rheumatoid arthritis. Rheumatoid Factor may appear transiently during various infections. The Rheumatoid Factor IgG ELISA is a diagnostic aid. A definite clinical diagnosis should not be based on the results of a single test, but should be made by the physician after all clinical and laboratory findings have been evaluated.